

Update on the Fast GC/Luminol System for NO₂ and PAN Measurements

Nancy A. Marley, Jeffrey S. Gaffney,
Environmental Research Division
Argonne National Laboratory

Abstract

The fast GC/Luminol instrument for the rapid simultaneous detection of NO₂ and PAN has been redesigned to be small and lightweight. It is contained in a standard 19 inch rack mountable case with a total weight of 35.5 pounds (not including laptop and carrier gas bottle). It is run from a laptop computer under a Labview environment. The software synchronizes injection and data collection, eliminating the need for hardware to actuate the injection valve. Both raw and integrated Data is stored automatically in spreadsheet format. The instrument can make use of Matheson Minimat gas cylinders (12in. X 2 in; 4 lbs) for both carrier and standard to minimize size and weight.

The reaction cell has been redesigned to minimize dead volume and the interface between the cell and PMT is constructed to minimize stray light. All DC power is supplied by one internal power supply.

The new design was field tested on the NOAA Twin Otter during the BRACE campaign (May, 2002) and operated unattended through most of the study. These results are compared to those obtained in central California during the CCOS field campaign (July, 2000) with the previous instrument design.


This new design allows for detection limits of 10 ppt for PAN with a PAN/NO₂ sensitivity Ratio of 1.4 depending on the Luminol conditions.

Analytical Determination of Reactive Nitrogen Species:


Chemiluminescent Reaction with Ozone:



If $\text{O}_3 \gg \text{NO}$ = NO detection



(also PANs, RNO_3 , NH_3 , HNO_3 ,
 N_2O_5 , ClNO_x , HONO)



(also HONO, NO_3 , HO_2NO_2)


PEROXYACYL NITRATES

- Organic oxidants having the general chemical structure $\text{RC}=\text{OO}-\text{O}-\text{NO}_2$
- Most common PANs are peroxyacetyl nitrate (PAN; $\text{R}=\text{CH}_3-$), peroxypropionyl nitrate (PPN; $\text{R}=\text{CH}_3\text{CH}_2-$), and peroxybutyryl nitrate (PBN; $\text{R}=\text{CH}_3\text{CH}_2\text{CH}_2-$).
- PANs are in thermal equilibrium with the peroxyacetyl radical ($\text{RC}=\text{O}-\text{OO}\cdot$) and $\text{NO}_2 \Rightarrow$ trapped peroxy radicals.
 - \Rightarrow indicator species of the photochemical age of an air parcel.
 - \Rightarrow important as a vehicle for long-range transport of NO_2 .
 - \Rightarrow leading to the formation of regional ozone and other oxidants.

Analytical Determination of Peroxyacyl Nitrates (PANs):




Gas Chromatography with Electron Capture Detection:



Also O₂, CFC's, H₂O have strong signals. – background limits analysis speed (typically 15 - 30 min.) before signal returns to baseline values.

Chemiluminescent Reaction with Luminol:




R = NO₂, PANs, H₂O₂, amino PAHs, metal ions, enzymes, carbohydrates

R can be selected by changing the pH, ionic strength or including a catalyst in the luminol solution.


Light intensity \propto R concentration

Chemiluminescent Emission Regions




Luminol Chemiluminescence is more sensitive than ozone chemiluminescence for reactive nitrogen species because the available photo multiplier detectors are more sensitive in the region of emission.

Fast GC with Luminol Detection




Hamamatsu HC135-1 Photon Counting Module

| | Specifications | Luminol Chemiluminescence |
|----------------------------|---------------------------------------|--|
| Spectral Range = | 300-650 nm | 375-550 nm |
| Peak Wavelength = | 400 nm | 425 nm |
| Counting Efficiency = | 22 % | |
| Dynamic Range = | 2×10^6 | |
| Linearity = | $0 \rightarrow 2 \times 10^6$ | Range = DL \rightarrow 160 ppb |
| Spectral Sensitivity = | 4.3×10^5 cps/pW | |
| Maximum Light Input = | 70 pW | $\cong 30$ Mcps |
| Dark Count = | 100 s^{-1} | $\cong 8$ ppt |
| Baseline Stability = | 11 %/ $^{\circ}\text{C}$ | $\Rightarrow 5^{\circ}\text{C} = 55$ cps drift |
| Response = | 0.1 %/ $^{\circ}\text{C}$ | |
| Optimum Integration Time = | 1 s | |
| Smallest Interval = | 10 ms | |
| Dimensions = | $1\frac{1}{4} \times 4\frac{3}{4}$ in | |
| Window Size = | 1 in | |



Labview Virtual Instrument Front Panel



Labview Virtual Instrument Block Diagram

Block Diagram

Set up Files

Set Inj. time

PMT
On/Off

Set Int. time


Stores both raw & integrated
Data to separate files.

Start Inj cycle


Receive data

Integrate peaks


Store to file




PAN & NO₂ Standards



NO₂ & PAN in Room Air



NO_2 , PAN, PPN, PBN Standard Separation
Luminol Chemiluminescence Detection




NOAA Twin Otter Inside Cabin

Laptop with Labview Software


PAN/NO₂ Instrument
Luminol Chemiluminescence



Calibration of the fast GC-luminol instrument for NO₂
at 1-sec, 0.5-sec, and 0.25-sec integration times.



Thermal decomposition of PAN samples to NO_2 , followed by reformation of PAN, as measured by the fast GC-luminol instrument.




The ratio of peak areas for PAN/ NO_2 gives a relative sensitivity of 1.4

Comparison of raw NO₂ Signals

Luminol Chemiluminescence - ■


Ozone Chemiluminescence with Catalytic Conversion - ▲

Ozone Chemiluminescence with Photolytic Conversion - ◆




Typical NO_2 and PAN profiles obtained over the Tampa area during the Bay Region Atmospheric Chemistry Experiment (BRACE)

(May 3, 2002).



Typical NO_2 and PAN profiles obtained over the Western Boundary during the Bay Region Atmospheric Chemistry Experiment (BRACE)

(May 15, 2002).



BRACE 05/02 Summary of Twin Otter PAN results.

| <u>Date</u> | <u>Max (ppt)</u> | <u>Avg (ppt)</u> | <u>Median (ppt)</u> | <u>Detection Limits</u> | <u>Objective</u> |
|-------------|------------------|------------------|---------------------|-------------------------|--------------------------|
| May 1 | 185 | 51 | 39 | 10 | Test |
| May 2 | 197 | 46 | 34 | 10 | Test |
| May 3 | 275 | 103 | 95 | 10 | Bay Survey |
| May 6 | 229 | 62 | 52 | 10 | Sea Breeze |
| May 8 | 270 | 51 | 46 | 10 | Sea Breeze |
| May 10 | 196 | 70 | 69 | 25 | Urban Plume |
| May 12 | 167 | 43 | 48 | 35 | Power Plant Plume |
| May 13 | 94 | 12 | <30 | 30 | Urban Plume |
| May 15 | 110 | 33 | 32 | 15 | Gulf |
| May 17 | 160 | 22 | <30 | 30 | Gulf/Ground Sites |
| May 20 | 188 | 54 | 50 | 20 | Urban Plume |
| May 22 | 129 | 40 | 35 | 15 | Eastern/Western Boundary |
| May 28 | 83 | 31 | 31 | 15 | Eastern Boundary |
| May 29 | 122 | 36 | 35 | 15 | Urban Plume |
| May 30 | 64 | 20 | 19 | 15 | Urban Plume |
| May 31 | 104 | 8 | <20 | 20 | Urban plume |
| June 1 | 176 | 61 | 65 | 50 | Power Plant Plume |

Detection Limits reported in ppt based on a concentration which gives a signal equal to 3 times the standard deviation of baseline values.
PMT replaced on May 21.

High detection limits on June 1 probably due to Luminol decay.

BRACE 05/02 Summary of Twin Otter NO₂ results.


| <u>Date</u> | <u>Max (ppb)</u> | <u>Avg (ppb)</u> | <u>Median (ppb)</u> | <u>Detection Limits</u> | <u>Objective</u> |
|-------------|------------------|------------------|---------------------|-------------------------|--------------------------|
| May 1 | 4.57 | 0.50 | 0.11 | 0.010 | Test |
| May 2 | 0.273 | 0.08 | 0.07 | 0.010 | Test |
| May 3 | 16.42 | 0.91 | 0.35 | 0.010 | Bay Survey |
| May 6 | 8.10 | 0.38 | 0.11 | 0.010 | Sea Breeze |
| May 8 | 20.44 | 0.30 | 0.07 | 0.010 | Sea Breeze |
| May 10 | 2.63 | 0.14 | 0.11 | 0.020 | Urban Plume |
| May 12 | 21.69 | 0.45 | 0.12 | 0.025 | Power Plant Plume |
| May 13 | 6.11 | 0.37 | 0.10 | 0.020 | Urban Plume |
| May 15 | 8.62 | 0.28 | 0.12 | 0.010 | Gulf |
| May 17 | 21.49 | 0.27 | 0.07 | 0.020 | Gulf/Ground Sites |
| May 20 | 25.71 | 0.68 | 0.17 | 0.015 | Urban Plume |
| May 22 | 4.00 | 0.19 | 0.09 | 0.010 | Eastern/Western Boundary |
| May 28 | 6.20 | 0.13 | 0.08 | 0.010 | Eastern Boundary |
| May 29 | 19.05 | 0.65 | 0.09 | 0.010 | Urban Plume |
| May 30 | 13.82 | 0.23 | 0.05 | 0.010 | Urban Plume |
| May 31 | 16.05 | 0.51 | 0.07 | 0.015 | Urban plume |
| June 1 | 17.11 | 0.58 | 0.20 | 0.036 | Power Plant Plume |

Detection Limits reported in ppt based on a concentration which gives a signal equal to 3 times the standard deviation of baseline values.


PMT replaced on May 21.

High detection limits on June 1 probably due to Luminol decay.

Typical NO₂ and PAN profiles obtained over the Fresno area during the Central California Oxidant Study (CCOS)
(July 9, 2000).



Typical NO_2 and PAN profiles obtained over the Western Boundary during the Central California Oxidant Study (CCOS)
(July 8, 2000).



CCOS 07/00 Summary of G1 PAN results.

| <u>Date</u> | <u>Max (ppt)</u> | <u>Avg (ppt)</u> | <u>Median (ppt)</u> | <u>Detection Limits</u> | <u>Objective</u> |
|-------------|------------------|------------------|---------------------|-------------------------|-------------------------|
| June 28 | 225 | 76 | <70 | 70 | Northern Boundary |
| July 5 | 1206 | 150 | 96 | 70 | Fresno Area, Low Passes |
| July 5 | 776 | 133 | 106 | 70 | Fresno Area |
| July 8 | 2969 | 71 | <70 | 70 | Western Boundary, AM |
| July 8 | 2908 | 91 | <70 | 70 | Western Boundary, PM |
| July 9 | 798 | 163 | 100 | 70 | Fresno Area, Low Passes |
| July 9 | 933 | 169 | 155 | 70 | Fresno Area |
| July 10 | 993 | 134 | 94 | 70 | Fresno Area, Low Passes |
| July 10 | 2965 | 231 | 211 | 70 | Fresno Area |
| July 11 | 3065 | 261 | 140 | 70 | Fresno Area, Low Passes |
| July 11 | 2157 | 308 | 191 | 70 | Fresno Area |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Detection Limits reported in ppt based on a concentration which gives a signal equal to 3 times the standard deviation of baseline values.

High detection limits due to use of air as carrier.

ACKNOWLEDGMENTS

This effort was supported by the U. S. Department of Energy (USDOE), Office of Science, Office of Biological and Environmental Research, Atmospheric Chemistry Program, under contract W-31-109-Eng-38.

We thank Mr. Peter Lunn (USDOE) for his continuing encouragement.

The aircraft data taken during the BRACE field campaign.
Was supported by the NOAA Air Resources Laboratory and the Florida Department of Environmental Protection.

CCOS 07/00 Summary of G1 NO₂ results.

| <u>Date</u> | <u>Max (ppb)</u> | <u>Avg (ppb)</u> | <u>Median (ppb)</u> | <u>Detection Limits</u> | <u>Objective</u> |
|-------------|------------------|------------------|---------------------|-------------------------|-------------------------|
| June 28 | 2.03 | 0.27 | 0.016 | 0.05 | Northern Boundary |
| July 5 | 49.80 | 1.60 | 0.12 | 0.05 | Fresno Area, Low Passes |
| July 5 | 19.00 | 0.77 | 0.26 | 0.05 | Fresno Area |
| July 8 | 52.73 | 0.97 | <0.05 | 0.05 | Western Boundary, AM |
| July 8 | 23.57 | 0.48 | 0.05 | 0.05 | Western Boundary, PM |
| July 9 | 28.91 | 1.55 | 0.14 | 0.05 | Fresno Area, Low Passes |
| July 9 | 17.94 | 0.54 | 0.29 | 0.05 | Fresno Area |
| July 10 | 15.34 | 0.65 | 0.12 | 0.05 | Fresno Area, Low Passes |
| July 10 | 9.89 | 0.68 | 0.41 | 0.05 | Fresno Area |
| July 11 | 53.59 | 3.53 | 0.195 | 0.05 | Fresno Area, Low Passes |
| July 11 | 35.84 | 1.77 | 0.343 | 0.05 | Fresno Area |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Detection Limits reported in ppt based on a concentration which gives a signal equal to 3 times the standard deviation of baseline values.

High detection limits due to use of air as carrier.